

Model 2781

THE FREEDOM BRIDGE

Free-Standing Rail System- No Wall or Ceiling Modifications



SUREHANDS®




- Perfect option for apartment or condo.
- Adjusts to fit the space available.
- Easy to install and move, can be a temporary or permanent solution.
- Can be flush mounted against the wall.
- Compatible with a fixed or portable motor.



Specialists in innovative and reliable lift & care systems

www.surehands.com

THE FREEDOM BRIDGE



● Technical specifications

GENERAL

Safe load: 529 lbs / 240 kg
 Rail: reinforced, aluminum white lacquered
 Side supports: aluminum, gray lacquered

DIMENSIONS

A Length:.....max 16.4 feet / 5000 mm (choose the rail length based on the available area)
B Height: max 8.2 feet/ 2500 mm (can be shortened based on the available height)
C Length foot: 1000 mm
D Height foot:..... 30 mm



Local Representative

Lift & Care® is a registered trademark of
 T.F. Herceg, Inc. dba SUREHANDS LIFT & CARE SYSTEMS.

© 21-11-2016 • Handi-Move • v1.1