




HANDI-MOVE MODEL 3200 + 3240

WHEELCHAIR-TO-WATER™ + SOCKET A



WHEELCHAIR-TO-WATER™



Dimensions (in inch/mm)

A	72.44 / 1840.00
B	51.38 / 1305.00
C	21.26 / 540.00
D	13.78 / 350.00
E	66.93 / 1700.00
F	21.65 / 550.00
G	19.69 / 500.00

Safe working load	300 Lbs / 135.00 kg
Batteries	2 x 12 Vdc; 6,5 - 7,2 Ah
Charger	separate

Diameter of castors	3.94 inch / 100.00 mm
Construction	stainless steel 316
Hand control	electrical IP67

SOCKET A




BODY SUPPORT SYSTEMS



Body Support

S	ref. 1025
S electrolytically coated (anticorrosive)	ref. 1027



+ Slipproof thigh supports ref. 1080



Standard Sling with moon shaped head support and padded thigh flaps

XS	ref. 3433
S	ref. 3430
M	ref. 3431
L	ref. 3432
XL	ref. 3434



+ Spreaderbar
ref. 2110 + 2111



Seat Sling (PVC)

One size	ref. 3473
----------	-----------

+ Spreaderbar ref. 2110



Bathing Sling (PVC)

XS	ref. 3476
S	ref. 3463
M	ref. 3464
L	ref. 3465
XL	ref. 3477

+ Spreaderbar ref. 2110

