
Ceiling Track Rails

Discreet and Efficient Transferring Solution



SUREHANDS®



- Customized track designs can provide single or multiple transfer locations. Combine with curves, switches and traverse track to accomplish transfers where needed.
- Rails can be installed on the ceiling or from wall to wall.
- Aesthetic design with warm white color.


Specialists in innovative and reliable lift & care systems

www.surehands.com


Ceiling Track Rails



● Dimensions (inch / mm)



Wall-to-wall rails



● Technical specifications

Safe working load..... 529 - 727.5 lbs / 240 - 330 kg
 Material..... Aluminum painted

Local Representative

Lift & Care® is a registered trademark of
 T.F. Herceg, Inc. dba SUREHANDS LIFT & CARE SYSTEMS.

© 10-04-2020 • Handi-Move • v1.0

SUREHANDS LIFT & CARE SYSTEMS
 982 County Route 1, Pine Island, NY 10969
 Tel: (800) 724-5305 • Fax: (845) 258-6634 • Website: www.surehands.com