
Traverse Rail

Full Room Coverage



SUREHANDS®




- Transfer anywhere in the room.
- Parallel rails can be mounted to the ceiling, or wall to wall.
- A manual version is easy to use with little effort from the caregiver.
- An electric version is available so the entire system can be used independently.
- Accessing an adjoining room is possible with the optional traverse connection.


Specialists in innovative and reliable lift & care systems

www.surehands.com

Traverse rail



● Dimensions (inch / mm)



● Technical specifications

Safe working load529 lbs / 240 kg
 MaterialAluminum painted, ELVZ steel, ABS

Local Representative

Lift & Care® is a registered trademark of T.F. Herceg, Inc. dba SUREHANDS LIFT & CARE SYSTEMS.

© 10-04-2020 • Handi-Move • v1.0