

Model 1030-1035

SureHands® Body Support



SUREHANDS®

The body support enables you to live more independently



- Lifts and transfers quickly, easily and safely
- Can be used independently or with an attendant
- Move from wheelchair to bed, bath, shower, toilet...
- Ideal for aquatic, hippo and ambulation therapies
- Leaves the seat entirely free - convenient for toilet use
- Comfortable and ergonomic sitting position
- Patented thigh supports designed for extra comfort
- Hygienic. Easy to clean

Specialists in innovative and reliable lift & care systems

www.surehands.com




SureHands® Body Support

● Dimensions (inch / mm)

	M	S	XS
A	31.1 / 790	27.95 / 710	26.38 / 670
B	18.9 / 480		


■ Size determined by lifting height requirements.



● Thigh supports

The thigh supports are key to the function of the Body Support while ensuring the user is in a comfortable and proper position.

Under the weight of the user, the self-locking feature holds the thigh supports in position. Once the weight is gone, they can easily be moved.



● Technical specifications

Material "arms"	316 stainless steel
Material "hands"	polyurethane
Maximum load	400 lbs

● Versions (M - S - XS)

1030 - Standard


1035 - For the swimming pool: electropolished

● Accessories


Utility band




Cushions for curve-around cups



Thigh supports in different dimensions and materials



Fabric covers for thigh supports



Local Representative

Lift & Care® is a registered trademark of T.F. Herceg, Inc. dba SUREHANDS LIFT & CARE SYSTEMS.

© 13-09-2018 • Handi-Move • v2.1