## Thigh supports



The key to making the body support self-adjust for you







- SureHands® offers a variety of thigh supports making it possible for many people to be able to use the SureHands® body support.
- Thigh supports fit securely under the thigh and distribute a person's body weight and pressure uniformly. They also ensure that the body maintains posture and is in balance during lifting.
- Thigh supports are made of a maintenance friendly rust-resistant material.
- The patented thigh support fix<sup>™</sup> designed for extra comfort.




## Thigh supports

## Versions and dimensions



Art. n°	Description	Color	A [inch/ mm]	B [inch/ mm]
10745	Standard slipproof small	black	6.30 / 160	5.12 / 130
1074M	Standard slipproof medium	black	8.66 / 220	7.09 / 180
1074L	Standard slipproof large	black	11.22 / 285	9.84 / 250
1074XL	Standard slipproof extra-large	black	12.2 / 310	9.84 / 250
1075M	Contour slipproof medium	black	8.66 / 220	5.91 / 150
10765	Muscular dystrophy slipproof small	black	6.30 / 160	5.12 / 130
1076M	Muscular dystrophy slipproof medium	black	8.66 / 220	7.09 / 180
10775	Standard small	grey	6.3 / 160	5.12 / 130
1077M	Standard medium	grey	8.66 / 220	7.09 / 180
1077XL	Standard extra-large	grey	12.2 / 310	9.84 / 250
1078M	Contour medium	grey	8.66 / 220	5.91 / 150
10795	Muscular dystrophy standard small	grey	6.3 / 160	5.12 / 130
1079M	Muscular dystrophy standard medium	grey	8.66 / 220	7.09 / 180



Local Representative

Lift & Care® is a registered trademark of T.F. Herceg, Inc. dba SureHands Lift & Care Systems.

© 04-09-2017 • Handi-Move • v3.0