

Model 2800-2805

2800 series motor track systems

Efficient and fast



SUREHANDS®



- The fastest track motor
- Suitable for intensive use
- Featherlight sideways movement
- To be used with slings and/or the SureHands® Body Support
- Motor and rail systems make it possible to access all areas
- Wireless control (optional)
- Back-up control on motor
- Emergency down and stop within easy reach

Lift & Care systems for people with reduced mobility

www.surehands.com




2800 series motor track systems

● Technical specifications

Safe working load	440 lbs / 200 kg - 660 lbs / 300 kg (118.11 inch / 3000 mm motorband: 352.74 lbs / 160 kg)
Batteries	2 x 12 V 9.5 Ah/ rechargeable
Standard motor band length	78.7 inch / 2000 mm (option: 118.11 inch / 3000 mm)
Rail System	aluminium, lacquered in white
Smooth stop and go	electronic
Cover	ABS
Emergency down	mechanical
Handcontrol 2800	watertight IP67 electric
Handcontrol 2805	infra-red
Number of lifts with fully charged batteries	160
Load and battery indications	visual and auditive
Charger	input 0.9 A 100-240 Vac,50-60 Hz Class II
.....	output 24 Vdc Max 1.3 A (9940) Max 1.2 A (2541)

● Examples of routes




● Versions*

2800
Lifting motor (up/down)
With electric handcontrol

2805 IR
Lifting motor (up/down)
With infra-red remote control

* also available with lateral motor (2810 L - 2815 L IR)

Local Representative



Lift & Care® is a registered trademark of
T.F. Herceg, Inc. dba SUREHANDS LIFT & CARE SYSTEMS.

© 30-11-2015 • Handi-Move • v1.3