

Model 2800-2805

Ceiling Lift

Ideal for Use with Caregiver



SUREHANDS®



- The fastest ceiling lift
- Lifts from floor
- Customized track systems make it possible to access all areas
- Accommodates SureHands® Body Support and Handi-Slings
- Fast up and down for efficient transfers
- Mechanical emergency down feature
- Infrared wireless remote control available
- Smooth lateral movement
- Automatic stop at any desired spot

Specialists in innovative and reliable lift & care systems

www.surehands.com

● Technical specifications

Safe working load	440 lbs / 200 kg (118.11 inch / 3000 mm motorband: 352.74 lbs / 160 kg)
Batteries	2 x 12 V 8 Ah / rechargeable
Standard motor band length	78.7 inch / 2000 mm (option: 118.11 inch / 3000 mm)
Rail System	aluminium, lacquered in white
Smooth stop and go	electronic
Cover	ABS
Emergency down	mechanical
Handcontrol 2800	watertight IP67 electric
Handcontrol 2805	infra-red
Number of lifts with fully charged batteries	160
Load and battery indications	visual and auditive
Charger	input 0.9 A 100-240 Vac, 50-60 Hz Class II
.....	output 24 Vdc Max 1.3 A (9940) Max 1.2 A (2541)

● Speeds


Speed raising

Load 5 kg	2.36 inch / 6.00cm/sec
Load 100 kg	2.28 inch / 5.80 cm/sec
Load 200 kg	1.81 inch / 4.60 cm/sec

Speed lowering

Load 5 kg	2.36 inch / 6.00 cm/sec
Load 100 kg	2.44 inch / 6.20 cm/sec
Load 200 kg	2.59 inch / 6.60 cm/sec

● Examples of routes



● Versions*

2800

Lifting motor (up/down)
With electric handcontrol

2805 IR

Lifting motor (up/down)
With infrared remote control

* also available with lateral motor (2810 L - 2815 L IR)

Lift & Care® is a registered trademark of
T.F. Herceg, Inc. dba SUREHANDS LIFT & CARE SYSTEMS.

© 21-01-2025 • Handi-Move • v.1.9

Local Representative