

Model 2840

# Portable Ceiling Motor

Multiple transfer locations with one motor



SUREHANDS®





- Install rail in rooms where transfers are needed (bed, bath, toilet, shower...).
- No doorway modifications needed.
- Light, compact, only 17 ½ lbs, yet powerful lifts 400 lbs.
- Economical: multiple transfers with one motor.
- Streamlined design fits home aesthetics.
- Economical: multiple transfers with one motor.
- Use with Body Support or slings.
- 10-year design life that you can count on.
- Option: transport trolley for the portable motor.

Specialists in innovative and reliable lift & care systems


[www.surehands.com](http://www.surehands.com)

# Portable Ceiling Motor

## Dimensions



A	18.89 inch - 480 mm
B	7.8 inch - 200 mm
C	8.27 inch - 210 mm
D	61.02 inch - 1550 mm




## Technical specifications

Weight .....	17.5 Lbs - 8 kg
Max. Lifting Capacity .....	400 Lbs - 182 kg
Batteries .....	2 x 12 VDC/3.2 Ah.
Charger: .....	input 220 V AC / output 24 V DC
Number of lifts with full batteries.....	50 x with 220 Lbs - 100 kg over a distance of 19.5 inch - 50 cm
Motor cover .....	ABS
Emergency down .....	Electric
Overhead rail composition.....	Aluminum alloy, white laquered

## With Handi-slings or SureHands® Body Support

Thanks to the integrated hooks on both sides of the motor and the swivel hook placed centrally at the bottom, the motor can be used with all Handi-slings or with the Body Support.



Option: Transport trolley

Local Representative

Lift & Care® is a registered trademark of  
T.F. Herceg, Inc. dba SUREHANDS LIFT & CARE SYSTEMS.

© 13-09-2023 • Handi-Move • v1.1

SUREHANDS LIFT & CARE SYSTEMS  
982 County Route 1, Pine Island, NY 10969  
Tel: (800) 724-5305 • Fax: (845) 258-6634 • Website: [www.surehands.com](http://www.surehands.com)