

Model 2845

Portable ceiling lift with “tarzan” system

“Tarzan” or “swing” under door between adjoining rooms



SUREHANDS®




- Access adjoining rooms without modifying doorways.
- Light, compact, only 17 ½ lbs, yet powerful: lifts to 400 lbs.
- Streamlined design fits home aesthetics.
- Economical: multiple transfers with one ceiling lift.
- Use with Body Support or slings.
- 10-year design life that you can count on.
- Option: transport trolley for the Portable ceiling lift.

Specialists in innovative and reliable lift & care systems

www.surehands.com

Portable ceiling lift with "tarzan" system



Dimensions



A	18.89 inch - 480 mm
B	7.8 inch - 200 mm
C	8.27 inch - 210 mm
D	61.02 inch - 1550 mm

Specifications

Weight	17.5 Lbs - 8 kg
Max. Lifting Capacity	400 Lbs - 182 kg
Batteries	2 x 12 VDC/3.2 - 3.6 Ah
Charger:	input 220 V AC / output 24 V DC
Number of lifts with full batteries.....	50 x with 220 Lbs - 100 kg over a distance of 19.5 inch - 50 cm
Motor cover	ABS
Emergency down	Electric
Overhead rail composition.....	Aluminum alloy, white laquered



With slings or Body Support

Thanks to the integrated hooks on both sides of the motor and the swivel hook placed centrally at the bottom, the Portable ceiling lift can be used with all slings or with the Handi-Move Body Support.

Lift & Care® is a registered trademark of
T.F. Herceg, Inc. dba SUREHANDS LIFT & CARE SYSTEMS.

© 28-04-2025 • Handi-Move • v1.3

Local Representative

