

Model 3370

Shower bath

Shower, bath and treatment in one



SUREHANDS®



- Unique concept: shallow bath, stretcher shower and nursing table.
- Electrically-adjustable sides, deep/shallow and height adjustable.
- Very comfortable and easy to maintain.
- Moves on casters, mobile lift can pass underneath.
- Removable side bars available.



Lift & Care systems for people with reduced mobility

www.surehands.com

● Dimensions and specifications

Maximum Safe Working Load 440.9 Lbs / 200 kg
 (including patient, accessories and water)
 Trolley mass 220.5 Lbs / 100 kgs

Trolley dimensions:

Trolley overall length 85.24 Inch / 2165 mm
 Trolley overall width 30.51 Inch / 775 mm
 Tub platform length 71.85 Inch / 1825 mm
 Tub platform width 18.90 Inch / 480 mm
 Wheel base 45.87 Inch / 1165 mm
 Wheel track 27.17 Inch / 690 mm

Adjustment range:

Maximum overall height 46.61 Inch / 1184 mm
 Trolley maximum height 43.66 Inch / 1109 mm
 Minimum overall height 24.76 Inch / 629 mm
 Trolley minimum height 21.81 Inch / 554 mm
 Maximum patient platform height 36.57 Inch / 929 mm
 Maximum tub depth 8.27 Inch / 210 mm






Electrical specifications:

Rechargeable battery 24V DC - 2.9 Ah
 Battery charger input 110 V AC
 Battery charge time: 6 hours (max)
 Number of full load cycles 30


Time to raise or lower patient platform 30 sec
 Time to raise or lower tub sides 30 sec

Accessories:

- Spare battery
- Foam wedge
- Option:
- Side bars in stainless steel



Local Representative



Lift & Care® is a registered trademark of
 T.F. Herceg, Inc. dba SUREHANDS LIFT & CARE SYSTEMS.

© 047-09-2015 • Handi-Move • v2.1