

Model 3370

Shower bath

Shower, bath and treatment in one



SUREHANDS®



- Unique concept: shallow bath, stretcher shower and nursing table.
- Electrically-adjustable sides, deep/shallow and height adjustable.
- Very comfortable and easy to maintain.
- Moves on casters, mobile lift can pass underneath.
- Removable side bars available.



Specialists in innovative and reliable lift & care systems

www.surehands.com

● Dimensions and specifications

Maximum Safe Working Load 440.9 Lbs / 200 kg
(including patient, accessories and water)

Trolley mass 220.5 Lbs / 100 kgs

Trolley dimensions:

Trolley overall length 85.24 Inch / 2165 mm

Trolley overall width 30.51 Inch / 775 mm

Tub platform length 71.85 Inch / 1825 mm

Tub platform width 18.90 Inch / 480 mm

Wheel base 45.87 Inch / 1165 mm

Wheel track 27.17 Inch / 690 mm

Adjustment range:

Maximum overall height 46.61 Inch / 1184 mm

Trolley maximum height 43.66 Inch / 1109 mm

Minimum overall height 24.76 Inch / 629 mm

Trolley minimum height 21.81 Inch / 554 mm

Maximum patient platform height 36.57 Inch / 929 mm

Maximum tub depth 8.27 Inch / 210 mm

Electrical specifications:

Rechargeable battery 24V DC - 2.9 Ah

Battery charger input 110 V AC

Battery charge time: 6 hours (max)

Number of full load cycles 30

Time to raise or lower patient platform 30 sec

Time to raise or lower tub sides 30 sec






Accessories:

Spare battery


Foam wedge

Option:

Side bars in stainless steel



Local Representative



Lift & Care® is a registered trademark of
T.F. Herceg, Inc. dba SUREHANDS LIFT & CARE SYSTEMS.

© 19-05-2016 • Handi-Move • v2.1

SUREHANDS LIFT & CARE SYSTEMS

982 County Route 1, Pine Island, NY 10969

Tel: (800) 724-5305 • Fax: (845) 258-6634 • Website: www.surehands.com